



**ADDENDUM TO
CHAMPION MANUFACTURING INC.
OPERATOR MANUALS**

**INFORMATION REGARDING INTEGRATED
OR STAND-ALONE USB OPTIONS**

ISSUE DATE: APRIL 21, 2020

Addendum information is to be added to Operator Manuals and published at a later date.
Refer to the Champion website for the most recent operator manual revisions.

CONTENTS

ITEM	PAGE
Contents	1
In-Service Information	2
Special notes	3
Power requirements	4
Product Labels and Symbols	5
Operation and General Cleaning of USB Charge Port	6
Power connections	7
Service information and Parts Identification	8
Recliner Schematics- Exploded parts view	9
End of Manual Addendum	12

Please see Operators Manual for all other chair information including, but not limited to:

- Special Notes (extended information)
- Chair Capacity Information
- Power Requirements (extended information)
- Receiving New Chairs
- Disposal of Equipment and Accessories
- Technical Description
- Product Labels (extended information)
- In-Service Procedures
- Operating Instructions (extended information)
- General Maintenance and Care (extended information)
- Schematics (extended information)
- Warranty and Service Information (extended information)
- Chair Options and Accessories



In-Service Information

It is very important to Champion that our customers review all the materials in the in-service package they receive.

By reviewing the in-service materials provided, you increase the useful life of the products by helping to prevent their improper use, which may lead to premature damage and wear. You are also helping to validate your warranty, as warranties cover material and workmanship failures, not misuse or abuse. Use these in-service tools for new employees and as a "refresher" for current employees.

The following resources are available for you to assist your staff and patients in operating the chairs properly.

Website: www.championchair.com

Operators Manuals

Find Operational Videos at: www.youtube.com/user/ChampionChair1992

Field Sales Reps

In House Sales Reps

Customer Care Team

Network of Service Technicians

In House Service Team

You should receive an Operator's Manual with each shipment of chairs you receive from Champion. If you would like information on any of these additional resources, please contact Champion's Customer Care Team (previously referred to as Champion's customer service department) at 800-998-5018.


NOTE: The Operator Manual may not be all inclusive of chair functions, options and/ or accessories.


NOTE: The information contained in this document is subject to change without notice.


SPECIAL NOTES

WARNING/CAUTION These terms refer to actions that may result in injury to your occupant or staff member and/ or damage to your product. Damage caused by improper operation of your chair is not covered by your warranty.

 **WARNING:** DO NOT operate this product without reviewing the in-service manual that demonstrates proper operation. This information accompanied your shipment

 **WARNING:** The power chair must be plugged into a properly grounded outlet and must not be modified in any way. If the 3-blade plug does not fit your outlet, have a new outlet installed that meets your local codes by an electrician. DO NOT use a 3-prong to 2-prong plug adapter. Consult an electrician if you are unsure that the outlets are properly grounded.

 **WARNING:** Keep hands, all other body parts and equipment away from moving parts that can cause pinch points! Do not reach to the inside of the chair unless the chair is unplugged from the outlet and the optional battery is disconnected to avoid injury if the chair position is changed. Caregiver assisting an occupant should always locate and safely place the occupant's extremities away from possible pinch points.

 **CAUTION:** Operating the recliner mechanism when the chair back is restricted and cannot move up and down freely may damage your mechanism. Be certain that the chair back as well as the push handle DO NOT make contact with the walls, counters or other furniture or equipment.

NOTE: The information contained in this document is subject to change without notice.

POWER REQUIREMENTS

Standard Power Requirements for Chair:


Electrical Requirements:	AC Power Supply: 120 VAC, 60 Hertz	
Rated Current Input:	Max. 3.7 Amps	
Power Output:	12VDC (Heat/Massage with or without integrated USB), 24VDC (Power recline), 5VDC (Standalone USB)	
Ambient Temperature:	+41° to +104°F (+5° to +40°C)	
Duty Cycle (Power recline)	10%, Max. 2 min./18min.	
IEC Protection Class:	Class II (Heat/Massage)	Class I (Power Recline)
Battery Duration (approx. times):	Backup battery system to remain plugged in. Approx. 15 recline cycles of the chair from battery.	
Battery Recharge Time:	4 Hours, Storage 10 Hours every 6 months Running battery to low power will cause damage.	

PRODUCT LABELS

Serial Number:



Name plate label
(Only where required)



Cleaning / Trend:




Cleaning



Trend

Standalone USB:



NOTE: Labels are not shown to actual size.

GRAPHIC SYMBOLS

Protective Earth Ground:	
Alternating Current / Direct Current:	
Do not Oil or Grease	
Keep Dry	
Attention: Consider Accompanying Documents	
Double Insulation	

NOTE: Symbols may not be shown to actual size.

USB CHARGE PORT- OPERATION



CAUTION: DO NOT allow moisture to pool in or around the USB Charging Port.


The USB Charging Port is strictly a charging port and does not transfer data. This is applicable to both the integrated and stand-alone USB options. Charging time of device will vary depending on charging cord and device type.

Operation:

Open the rubber dust cover and insert USB charging cord.

The blue light on the port glows when the charging port is in operation.

Remove the charging cord and replace dust cover when port is not in use.



Note: Some components may have a different appearance than shown in the picture due to multiple models of chairs. Wires, fasteners, etc. may not be labeled or shown. The serial number of the chair informs Customer Care of the proper components for that chair.

USB CHARGE PORT- GENERAL CLEANING



CAUTION: DO NOT allow moisture to pool in or around the USB Charging Port.

The rubber dust cover does not provide a seal. Care should be taken while cleaning around the port as to not let moisture enter the electrical components.

If dust or debris is noticed the USB port may be cleaned using compressed air to blow the particles out.

POWER CONNECTIONS

Stand- alone USB


The stand-alone USB port is to be plugged into a properly grounded standard wall outlet. (Stand-alone USB outlet plug is not shown)

Integrated USB (into a power only chair)

The integrated USB port is to be plugged into a properly grounded standard wall outlet. (Integrated USB outlet plug is not shown)

Integrated USB (into a power, heat, and/ or massage chair)

The integrated USB port to a powered heat and/or massage chair is plugged into the Control box- Accessory Connector Plug. The Control box is then plugged into a properly grounded standard wall outlet by utilizing the AC Connector. (AC connector cord with plug is not shown)



Note: Some components may have a different appearance than shown in the picture due to multiple models of chairs. Wires, fasteners, etc. may not be labeled or shown. The serial number of the chair informs Customer Care of the proper components for that chair.

SERVICE INFORMATION

The mission of the Customer Care Team is to get your chair up and running as quickly as possible. It is critical that the Customer Care Team knows what product you have, and exactly what is wrong with the product. If you have questions or problems, you should never hesitate to call for assistance: 800-998-5018.

The most timely and cost-effective way for your chair to be repaired is for the Champion's Customer Care Team to work with your maintenance department or equipment technician.

Determining the Problem

What is wrong with the chair should be determined by troubleshooting. The Customer Care Team will assist you with this by asking you questions about the chair function.

Serial Number

The chair serial number identifies the precise configuration of your chair; this is critical to receiving correct components and instructions. **This number is required to process your request.**

The serial number is located in the back of the chair on the lower left side on the label entitled Champion Manufacturing -Serial #xxxxxx.

No parts will be shipped without a Serial Number to which they will be installed.

PARTS IDENTIFICATION

To identify worn or damaged components please refer to appropriate product schematics.

To obtain repair part numbers refer to the parts listing key using the schematics page and item number.


Parts orders may be placed by using the convenient fax order form in this manual or by calling Champion's Customer Care 800-998-5018.

When placing an order by phone you will be asked to provide:

- your name and facility name
- your phone, fax number, email address
- the serial number of your product, and
- the nature of your problem

Having the above information available at the time you call will expedite the process. In order to provide prompt, accurate service it may be necessary to request further information about the chair function to accurately define the problem.

SCHEMATICS FOR FIXED ARM CHAIRS




FIXED ARM CHAIR KEY COMPONENT LIST

PG/PART	DESCRIPTION	PG/PART	DESCRIPTION
09-01	COMPLETE FRAME BASE ASSEMBLY	09-05	UPHOLSTERED ARM INSERT (RIGHT WITH ACTUATOR AND USB SHOWN)
09-02	TABLE BRACKET SET	09-06	USB PORT WITH COVER (STAND-ALONE SHOWN)
09-03	TABLE TOP	09-07	ARM ASSEMBLY – FULL (RIGHT ACTUATOR SHOWN)
09-04	RECLINE MECHANISM ACTUATOR HANDLE		

Note: Some components may have a different appearance than shown in the picture due to multiple models of chairs. Wires, fasteners, etc. may not be labeled or shown. The serial number of the chair informs Customer Care of the proper components for that chair.

SCHEMATICS FOR SWING ARM CHAIRS



SWING ARM CHAIR KEY COMPONENTS LIST (ON NEXT PAGE)

SWING ARM CHAIR KEY COMPONENTS LIST

PG/PART	DESCRIPTION	PG/PART	DESCRIPTION
10-01	MECHANISM WITH ALL UPHOLSTERED PARTS ATTACHED- FOR ILLUSTRATION	10-14	"C" HINGE ARM BRACKET
10-02	LEFT HAND SWING ARM ASSBMELY	10-15	3/8 X 16/9" HEX BOLT
10-03	SWING ARM FRAME BASE	10-16	"C" HINGE BASE BRACKET
10-04	TREAD FOR PULL OUT FOOT TRAY	10-17	RING COTTER
10-05	PULL OUT FOOT TRAY	10-18	CAP RELEASE HANDLE
10-06	CASTER - FRONT TOTAL LOCK	10-19	RELEASE HANDLE
10-07	CASTER - STEERING	10-20	RELEASE COVER
10-08	CASTER - SWIVEL	10-21	LATCH ASSEMBLY STRIKE – RIGHT
10-09	TRACK SET FOR PULL OUT FOOT TRAY	10-22	LATCH ASSEMBLY MOUNT – RIGHT
10-10	RELEASE ASSEMBLY	10-23	LATCH ASSEMBLY
10-11	TABLE TOP	10-24	LATCH ASSEMBLY MOUNT – LEFT (NOT SHOWN)
10-12	TABLE BRACKET SET	10-25	LATCH ASSEMBLY STRIKE – LEFT (NOT SHOWN)
10-13	RIGHT HAND ARM ASSEMBLY	10-26	USB PORT WITH COVER (STAND-ALONE SHOWN)

Note: Some components may have a different appearance than shown in the picture due to multiple models of chairs. Wires, fasteners, etc. may not be labeled or shown. The serial number of the chair informs Customer Care of the proper components for that chair.

END OF ADDENDUM